



Thermal quality control on domestic appliances



Motorcycle break testing



# FLIR A655sc

High resolution science grade LWIR camera

The A655sc helps engineers, researchers, and scientists see and accurately quantify thermal patterns, leakage, dissipation, and other heat-related factors in equipment, products, and processes in real time.

# **IMAGE QUALITY AND THERMAL SENSITIVITY**

FLIR A655sc is equipped with an uncooled, maintenance free, Vanadium Oxide (VoX) microbolometer detector that produces thermal images of 640 x 480 Pixels. These pixels generate crisp and clear detailed images that are easy to interpret with high accuracy. The FLIR A655sc will make temperature differences as small as 50 mK clearly visible.

## **HIGH SPEED WINDOWING**

The FLIR A655sc provides 14-bit data up to 50 frames per second at full frame  $640 \times 480$  resolution. It has a high speed windowing function that increases the output frame rate up to 200 Hz at a  $640 \times 120$  pixel window.

## GIGE VISION™ STANDARD COMPATIBILITY

GigE Vision allows fast image transfer using low cost standard cables up to 100 meters. With GigE Vision, hardware and software from different vendors can integrate seamlessly over gigabit ethernet connections.

# GENICAM™ PROTOCOL SUPPORT

GenlCam creates a common application programming interface (API) for cameras regardless of the interface technology or features implemented. Because the API for GenlCam cameras will always be the same, cameras like the A655sc camera can be easily integrated into third party software.

#### SOFTWARE

FLIR A655sc camera works seamlessly with FLIR ResearchIR Max software enabling intuitive viewing, recording and advanced processing of the thermal data provided by the camera.

## MATHWORKS® MATLAB

Control and capture data directly into MathWorksR Matlab software for advanced image analysis and processing.

## **KEY FEATURES**

- Uncooled microbolometer: 640 x 480 pixels
- Gigabit ethernet and usb interface
- \* Close-up and telephoto lenses available
- ResearchIR max software included
- \* Matlab compatible



# **Imaging Specifications**

| System Overview  | FLIR A655sc   |
|--|---|
| Detector Type  | Uncooled Microbolometer   |
| Spectral Range   | 7.5 – 14.0 µm   |
| Resolution   | 640 × 480   |
| Detector Pitch   | 17 μm   |
| NETD   | <30 mK  |
| Imaging  |   |
| Time Constant  | <8 ms   |
| Frame Rate (Full Window)   | 50 Hz   |
| Subwindow Mode   | User-Selectable 640 × 240 or 640 × 120  |
| Maximum Frame Rate<br>(@ Min. Window)  | 200 Hz (640 × 120)  |
| Dynamic Range  | 14-bit  |
| Digital Data Streaming   | Gigabit Ethernet (50/100/200 Hz) USB (25/50/100 Hz)   |
| Command and Control  | Gigabit Ethernet, USB   |
| Measurement  |   |
| Standard<br>Temperature Range  | -40°C to 150°C (-40°F to 302°F)<br>100°C to 650°C (212°F to 1,202°F)  |
| Optional<br>Temperature Range  | Up to 2,000°C (3,632°F)   |
| Accuracy   | ±2°C or ±2% of Reading  |
| Optics   |   |
| Camera f/#   | f/1.0   |
| Available Lenses   | 6.5 mm (80°), 13.1 mm (45°), 24.6 mm (25°), 41.3 mm<br>(15°), 88.9 mm (7°)  |
| Focus  | (13 ), 66.3 11111 (7 )  |
|  | Automatic or Manual (Motorized)   |
| Close-up / Microscopes   | ( ) ( )   |
| Close-up / Microscopes Image Presentation  | Automatic or Manual (Motorized)   |
|  | Automatic or Manual (Motorized)   |
| Image Presentation   | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm   |
| Image Presentation  Digital Data   | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm   |
| Image Presentation Digital Data General Operating  | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm Via PC Using ResearchIR Software  |
| Image Presentation Digital Data General Operating Temperature Range Storage  | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm Via PC Using ResearchIR Software  -15°C to +50°C (+5°F to 122°F)  |
| Image Presentation Digital Data General Operating Temperature Range Storage Temperature Range                                      | Automatic or Manual (Motorized)  Close-up 25 μm, 50 μm, 100 μm  Via PC Using ResearchIR Software  -15°C to +50°C (+5°F to 122°F)  -40°C to 70°C (-40°F to 158°F)  |
| Image Presentation Digital Data General Operating Temperature Range Storage Temperature Range Encapsulation                        | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm  Via PC Using ResearchIR Software  -15°C to +50°C (+5°F to 122°F)  -40°C to 70°C (-40°F to 158°F)  IP 30 (IEC 60529)  |
| Image Presentation Digital Data General Operating Temperature Range Storage Temperature Range Encapsulation Bump / Vibration       | Automatic or Manual (Motorized)  Close-up 25 μm, 50 μm, 100 μm  Via PC Using ResearchIR Software  -15°C to +50°C (+5°F to 122°F)  -40°C to 70°C (-40°F to 158°F)  IP 30 (IEC 60529)  25 g (IEC 60068-2-29) / 2 g (IEC 60068-2-6)                              |
| Image Presentation Digital Data General Operating Temperature Range Storage Temperature Range Encapsulation Bump / Vibration Power | Automatic or Manual (Motorized) Close-up 25 μm, 50 μm, 100 μm  Via PC Using ResearchIR Software  -15°C to +50°C (+5°F to 122°F)  -40°C to 70°C (-40°F to 158°F)  IP 30 (IEC 60529) 25 g (IEC 60068-2-29) / 2 g (IEC 60068-2-6)  12/24 VDC, 24 W Absolute Max. |

Power Connector,
Screw Terminal
2-pole: 10-30 VDC, <10W

Control and image streaming

USB H2 Connector:
Camera control and image streaming

Digital I/O Connector,

Screw Terminal 6-pole: Digital Out: 2 outputs, opto-isolated, 10-30V supply, 100mA. Digital In: 2 inputs, opto-isolated, 10-30 V.

#### PORTLAND

Corporate Headquarters FLIR Systems, Inc. 27700 SW Parkway Ave. Wilsonville, OR 97070 USA PH: +1 866.477.3687

### **BELGIUM**

FLIR Systems Trading Belgium BVBA Luxemburgstraat 2 2321 Meer Belgium PH: +32 (0) 3665 5100

#### SWEDEN

FLIR Systems AB Antennvägen 6, PO Box 7376 SE-187 66 Täby Sweden PH: +46 (0)8 753 25 00

www.flir.com NASDAQ: FLIR

# NASHUA

FLIR Systems, Inc. 9 Townsend West Nashua, NH 06063 USA PH: +1 603.324.7611

#### UK

FLIR Systems UK 2 Kings Hill Avenue Kings Hill West Malling - Kent ME19 4AQ United Kingdom PH: +44 (0)1732 220 011

Specifications are subject to change without notice @Copyright 2014, FLIR Systems, Inc. All other brand and product names are trademarks of their respective owners. The images displayed may not be representative of the actual resolution of the camera shown. Images for illustrative purposes only. [Created 08/14]

